

Phigenics *Legionella* Testing Program For Value Added Resellers (VAR)

Phigenics Analytical Services Laboratory

For questions, please contact Anne Burt, direct line: 651.200.4612

3S701 West Street, Warrenville, IL 60555

Main Office Phone: 844.850.4087

Tests and Test Combinations

Phigenics <i>Legionella</i> Test Products (w/ Part Number)	Test Method	Bacteria Type	Turnaround Time (once received at lab)
PVT TimeZero™ (Part No. MB-100)	TimeZero™ Culture Method	<i>Legionella pneumophila</i> serogroup 1 <i>Legionella pneumophila</i> serogroups 2-14 <i>Legionella</i> species Total Heterotrophic Aerobic Bacteria (THAB)	4 Business Days
PVT with ISO (Part No. MB-300)	TimeZero™ Culture Method	<i>Legionella pneumophila</i> serogroup 1 <i>Legionella pneumophila</i> serogroups 2-14 <i>Legionella</i> species Total Heterotrophic Aerobic Bacteria (THAB)	4 Business Days
	ISO 11731 Spread Plate Culture Method	<i>Legionella pneumophila</i> serogroup 1 <i>Legionella pneumophila</i> serogroups 2-14 <i>Legionella</i> species	10-14 Business Days
PVT PREMIUM (Part No. MB-210)	TimeZero™ Culture Method	<i>Legionella pneumophila</i> serogroup 1 <i>Legionella pneumophila</i> serogroups 2-14 <i>Legionella</i> species Total Heterotrophic Aerobic Bacteria (THAB)	4 Business Days
	ISO 11731 Spread Plate Culture Method	<i>Legionella pneumophila</i> serogroup 1 <i>Legionella pneumophila</i> serogroups 2-14 <i>Legionella</i> species	10-14 Business Days
	Next Day LegiPlex™ PCR Method (a molecular marker negative screen)	Real-time PCR detection of <i>Legionella pneumophila</i> sg 1, <i>Legionella pneumophila</i> sg 2-14, and <i>Legionella</i> genus	1 Business Day

FOR MORE INFORMATION, VISIT:
INFO.PHIGENICS.COM/VAR

Tests and Test Combinations (cont)

Phigenics <i>Legionella</i> Test Products (w/ Part Number)	Test Method	Bacteria Type	Turnaround time (once received at lab)
PVT SELECT (Part No. MB-600)	TimeZero™ Culture Method	<i>Legionella pneumophila</i> serogroup 1 <i>Legionella pneumophila</i> serogroups 2-14 <i>Legionella</i> species Total Heterotrophic Aerobic Bacteria (THAB)	4 Business Days
	Next Day LegiPlex™ PCR Method (a molecular marker negative screen)	Real-time PCR detection of <i>Legionella pneumophila</i> sg 1, <i>Legionella pneumophila</i> sg 2-14, and <i>Legionella</i> genus	1 Business Day
PVT Next Day LegiPlex™ PCR (Part No. MB-1430)	Next Day LegiPlex™ PCR Method (a molecular marker negative screen)	Real-time PCR detection of <i>Legionella pneumophila</i> sg 1, <i>Legionella pneumophila</i> sg 2-14, and <i>Legionella</i> genus	1 Business Day
PVT Standard ISO Spread Plate Legionella Test (Part No. MB-1010)	ISO 11731 Spread Plate Culture Method	<i>Legionella pneumophila</i> serogroup 1 <i>Legionella pneumophila</i> serogroups 2-14 <i>Legionella</i> species	10-14 Business Days
PVT Next Day LegiPlex™ PCR with Standard ISO Spread Plate Legionella Test (Part No. MB-1120)	Next Day LegiPlex™ PCR Method (a molecular marker negative screen)	Real-time PCR detection of <i>Legionella pneumophila</i> sg 1, <i>Legionella pneumophila</i> sg 2-14, and <i>Legionella</i> genus	1 Business Day
	ISO 11731 Spread Plate Culture Method	<i>Legionella pneumophila</i> serogroup 1 <i>Legionella pneumophila</i> serogroups 2-14 <i>Legionella</i> species	10-14 Business Days
PVT Swab Legionella Test (Part No. MB-1010)	<i>Legionella</i> Swab Culture Analysis (qualitative)	<i>Legionella pneumophila</i> serogroup 1 <i>Legionella pneumophila</i> serogroups 2-14 <i>Legionella</i> species	10-14 Business Days
PVT THAB Total Heterotrophic Aerobic Bacteria (Part No. MB-1010)	TimeZero™ Total Bacteria Culture Method	Total Heterotrophic Aerobic Bacteria	4 Business Days

The safety of your customers' water systems is one of their utmost concerns. When they are ready to validate their water systems, they want to be assured that their *Legionella* results are accurate and are received in a timely and efficient manner with comprehensive reporting. As an authorized Phigenics VAR, you have access to *Legionella* test services that are the most accurate, fastest and most comprehensive in the industry – allowing you to provide real value to your customers.